

ScleroSafe™



The beauty of treating varicose veins

A new technique that simplifies the procedure and shortens the treatment process

The ScleroSafe platform is an advanced non-thermal, non-tumescent technology solution that eliminates the need for surgery or costly capital equipment, while minimizing the pain and recovery time associated with conventional treatment.

ScleroSafe's unique, inverse action dual syringe (also known as the DPS) injects Asclera® into the vein through one syringe to induce chemical ablation while simultaneously aspirating blood through the second syringe – all within a single motion and with single-hand operation. The aspiration of the blood through the double lumen catheter results in a mechanically inducing negative pressure just behind the injected sclerosant agent.

ScleroSafe is used with the 5cc\10cc Dual Syringe, a specially designed 5F 150mm long or 5F 350mm long double lumen catheter and a dedicated guide wire.

For more information about ScleroSafe and how to order, contact us at asclera@methapharm.com or by phone at 1-833-766-8346 (VEIN).



INDICATION: ScleroSafe Peripheral Venous Aspiration & Infusion Kit with Dual Procedure Syringe (ScleroSafe) is intended for the delivery of Asclera®, an FDA-approved sclerosant, in the treatment of varicosities in superficial veins with a diameter of 2 to 3mm. ScleroSafe is a sterile, single use device designed to provide transient access to the peripheral venous system.

For complete ScleroSafe contraindications and warnings, please consult Package Instructions for Use. For Asclera indications, contraindications, adverse effects and warnings, please consult prescribing information, available at www.asclera.com.

Distributed by Methapharm, Inc. Copyright © 2024 Methapharm, Inc. All rights reserved. ScleroSafe is a trademark of VVT Medical.